



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Polyethylene Aprons

**Product Code :**  
LBNY-00017-11000482



**Description :**

High-quality virgin material combines flexibility and resistance to chemicals, fats and oils.  
Disposable, provide protection for light- and medium-duty applications  
Easy to clean  
Food preparation, processing and handling; food service/ catering; janitorial and maintenance work.  
Embossed  
Conveniently packaged  
600/Cs.

**Technical Specification :**

Polyethylene Aprons



**Laboratory Instrument India**