



Email : sales@laboratoryinstrumentindia.com

Product Name :
Polyethylene Aprons

Product Code :
LBNY-00017-11000482



Description :

High-quality virgin material combines flexibility and resistance to chemicals, fats and oils.
Disposable, provide protection for light- and medium-duty applications
Easy to clean
Food preparation, processing and handling; food service/ catering; janitorial and maintenance work.
Embossed
Conveniently packaged
600/Cs.

Technical Specification :

Polyethylene Aprons



Laboratory Instrument India