



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
105 White Coveralls

**Product Code :**  
LBNY-00017-11000489



**Description :**

Ideal for use in controlled environments to protect against particulates and other contaminants. Spinning process creates barrier to particles, micro-organisms, and non-hazardous liquids. Light and durable coveralls provide coverage, including feet, while staying comfortable. Use in controlled environments with high standards for particle and microbiological contamination control.

- Clean processed and sterile
- Serged seams
- Standard hood, elastic opening
- Elastic wrists and ankles
- Lightweight, comfortable and durable
- Lowest level particle shedding
- Recyclable
- Attached thumb loops
- Pharmaceutical, biotech, medical device, electronic settings.
- Zipper closure
- Attached boots, PVC sole.

**Technical Specification :**

105 White Coveralls

---

# Laboratory Instrument India

