



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Denim Shop Aprons, Two-Compartment Bib Pocket

**Product Code :**  
LBNY-00017-11000490



**Description :**

Non pre-laced polyester ties for neck and waist  
With metal grommets  
Available in choice of several pocket styles  
Preshrunk 100% cotton, 14 oz. denim minimizes shrinkage  
Color: Blue.  
Measure 36L x 28 in.W (91 x 71cm).

**Technical Specification :**

Denim Shop Aprons, Two-Compartment Bib Pocket



**Laboratory Instrument India**