



Email : sales@laboratoryinstrumentindia.com

Product Name :
Denim Shop Aprons, Two-Compartment Bib Pocket

Product Code :
LBNY-00017-11000490



Description :

Non pre-laced polyester ties for neck and waist
With metal grommets
Available in choice of several pocket styles
Preshrunk 100% cotton, 14 oz. denim minimizes shrinkage
Color: Blue.
Measure 36L x 28 in.W (91 x 71cm).

Technical Specification :

Denim Shop Aprons, Two-Compartment Bib Pocket



Laboratory Instrument India