



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Digital Stop Watch

**Product Code :**  
LBNY-0003-1010008



**Description :**

Digital Stop Watch

**Technical Specification :**

Type: Digital.

Time Display: Hour, Minute and Seconds.

Calendar Display: Day, Date & Month.Alarm.

Press Model Button once, if the stopwatch is on Time Mode, to bring it to Stopwatch mode.

Press Start / Stop Button to start the Stopwatch.

Press Start/Stop button again to stop the Stopwatch.



**Laboratory Instrument India**