



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

400 Comfort Fit Design Disposable Coverall with Respirator Fit Hood, Elastic Wrists and Ankles, NAFTA

**Product Code :**

LBNY-00017-11000498



**Description :**

400 coveralls with respirator fit hood, elastic wrists and ankles provide the protection and durability you need so you can focus on the task at hand. Excellent abrasion resistance due its unique, nonwoven material  
Can be used in a variety of applications including general maintenance/operations, asbestos and lead abatement, mold remediation, and environmental cleanups.  
400 provides light-weight inherent barrier protection against hazardous dry particles and aerosols, and nonhazardous light liquid splash.

**Technical Specification :**

400 Comfort Fit Design Disposable Coverall with Respirator Fit Hood, Elastic Wrists and Ankles, NAFTA



**Laboratory Instrument India**