



Email : sales@laboratoryinstrumentindia.com

Product Name :
Iso Clean 265 Series Lab Coats

Product Code :
LBNY-00017-11000500



Description :

Serged seams have multiple interlocking threads that are sewn around the raw edges of garment material to create a strong, stress-resistant seam.
Knit cuff for tighter fit at wrist.
Three front pockets.
A-line silhouette.
Front zipper closure for easy donning and doffing.
Iso Clean lab coats provide an inherent barrier to particles, microorganisms, and non-hazardous light liquid splash in controlled environments.
Raglan sleeve design for greater range of motion and comfort.
Bulk packaged in double transparent poly liners.

Technical Specification :

Iso Clean 265 Series Lab Coats



Laboratory Instrument India

