



Email : sales@laboratoryinstrumentindia.com

Product Name :

Iso Clean 265 Series Lab Coats

Product Code :

LBNY-00017-11000500



Description :

Serged seams have multiple interlocking threads that are sewn around the raw edges of garment material to create a strong, stress-resistant seam.

Knit cuff for tighter fit at wrist.

Three front pockets.

A-line silhouette.

Front zipper closure for easy donning and doffing.

Iso Clean lab coats provide an inherent barrier to particles, microorganisms, and non-hazardous light liquid splash in controlled environments.

Raglan sleeve design for greater range of motion and comfort.

Bulk packaged in double transparent poly liners.

Technical Specification :

Iso Clean 265 Series Lab Coats



Laboratory Instrument India

