



Email : sales@laboratoryinstrumentindia.com

Product Name :
68-4000 Green Coveralls with Hood

Product Code :
LBNY-00017-11000503



Description :

Protection: Permeation tested against over 190 chemicals, including chemical warfare agents; fully encapsulated versions also available

Comfort: Textile like interior improves wearer acceptance

Design: Elastic waist, wrists and ankles, double zip closure, double cuffs with knitted inner cuff, ultrasonically welded and taped seams, 2 piece respirator fit hood

Provide an exceptional barrier against organic and inorganic chemicals as well as biological agents.

Green Coveralls with Hood incorporate unique design features that enhance wearer comfort and safety.

Applications: Chemicals handling/transportation, oil-based mud protection, hazardous waste remediation, sewage purification installations, industrial/tank cleaning, emergency response (Level B), pharmaceutical, mining, agriculture.

Technical Specification :

68-4000 Green Coveralls with Hood



Laboratory Instrument India

