



Email : sales@laboratoryinstrumentindia.com

Product Name :
125 Coveralls

Product Code :
LBNY-00017-11000512



Description :

Useful when handling many wet and dry chemicals, and certain hazardous materials such as PCBs
Taped-Seamed models are sewn and taped to provide chemical resistance against heavy liquid splashes and tough seam stress
Laydown collar
Effective protection against a range of chemicals. Excellent for waste management, PCB clean-up, hazardous response, and nuclear environments.
Bound-Seamed models are tightly sewn with a reinforced outer binding to increase strength and particle barrier
Elastic wrists
Elastic ankles.
Front zipper closure
Storm flap with tape closure.

Technical Specification :

125 Coveralls



Laboratory Instrument India

