



Email : sales@laboratoryinstrumentindia.com

Product Name :
500 Coverall with Hood

Product Code :
LBNY-00017-11000514



Description :

500, model Coverall with hood. Serged external seams. Elastic wrists, ankles and face. Elastic waist (glued-in). Storm Flap. White.

500 coveralls offer enhanced protection against a range of threats, from water based chemical liquid splashes, asbestos, pharmaceutical and radioactive particulates to biological hazards.

Technical Specification :

500 Coverall with Hood



Laboratory Instrument India