



Email : sales@laboratoryinstrumentindia.com

Product Name :
182B Coveralls, Bulk

Product Code :
LBNY-00017-11000517



Description :

Bulk packaged in double transparent poly liners.
Made using a patented flash spinning process
Elastic openings for tighter fit at wrist and ankle
Bound neck for lower particle shedding
Front zipper closure for easy donning and doffing
Provides an inherent barrier to particles, microorganisms, and non-hazardous light liquid splash
Iso Clean delivers an ideal balance of protection, durability, and comfort. Bound seams are covered with garment fabric to reinforce the seam and to reduce the potential for particle penetration
Raglan sleeve design for greater range of motion and comfort.

Technical Specification :

182B Coveralls, Bulk



Laboratory Instrument India