



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
TK Training Suits

**Product Code :**  
LBNY-00017-11000518



**Description :**

Serged Seams. Lime Yellow. Clearly Labeled as a Training Suit. Size 10000 fully encapsulated training suits offer an economical and effective option for HazMat training. 10000 Encapsulated Training Suit. Expanded Back, Rear Entry. Extra-Wide Visor, 1 Layer: 20 mil PVC. Attached Butyl Gloves. Attached Socks with Outer Boot Flaps. Storm Flap. Same fabric and garment design as standard 10000 suits. Multiple interlocking threads are sewn around the raw edges of garment material to create a strong, stress-resistant seam

**Technical Specification :**

TK Training Suits



**Laboratory Instrument India**