



Email : sales@laboratoryinstrumentindia.com

Product Name :
3 Boot Covers (10 in. Tall)

Product Code :
LBNY-00017-11000519



Description :

Spunbonded polypropylene substrate coated with a low-density polyethylene film for enhanced skid resistance and barrier performance.

Spun bond polypropylene boot covers with elastic top provides coverage up to calf.

3 sole material provides enhanced skid-resistance

One size fits most.

Technical Specification :

3 Boot Covers (10 in. Tall)



Laboratory Instrument India