



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
3 Boot Covers (10 in. Tall)

**Product Code :**  
LBNY-00017-11000519



**Description :**

Spunbonded polypropylene substrate coated with a low-density polyethylene film for enhanced skid resistance and barrier performance.  
Spun bond polypropylene boot covers with elastic top provides coverage up to calf.  
3 sole material provides enhanced skid-resistance  
One size fits most.

**Technical Specification :**

3 Boot Covers (10 in. Tall)



**Laboratory Instrument India**