



Email : sales@laboratoryinstrumentindia.com

Product Name :
Safe Wear High Performance Lab Coat

Product Code :
LBNY-00017-11000520



Description :

High quality ultra-soft & light weight SMS fabric (Weight 45 g/m2 for Jacket & Lab Coat)
Designed to fit properly with knitted cuffs & collar
High-quality front snap closure buttons
Safe Wear premium protective apparel is made from top-quality materials and designed to provide reliable barrier protection without sacrificing comfort or style.
Fluid-resistant material meets AAMI Level 2
Available in a wide range of sizes and in fashionable colors
12 units in every pack.

Technical Specification :

Safe Wear High Performance Lab Coat



Laboratory Instrument India