



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Safe Wear High Performance Lab Coat

**Product Code :**  
LBNY-00017-11000521



**Description :**

High quality ultra-soft & light weight SMS fabric (Weight 45 g/m2 for Jacket & Lab Coat)  
Fluid-resistant material meets AAMI Level 2  
Designed to fit properly with knitted cuffs & collar  
High-quality front snap closure buttons  
Safe Wear premium protective apparel is made from top-quality materials and designed to provide reliable barrier protection without sacrificing comfort or style.  
Available in a wide range of sizes and in fashionable colors  
12 units in every pack.

**Technical Specification :**

Safe Wear High Performance Lab Coat



**Laboratory Instrument India**