



Email : sales@laboratoryinstrumentindia.com

Product Name :
11-280 Sleeves

Product Code :
LBNY-00017-11000522



Description :

INTERCEPT Cut Resistance Technology provides best-in-class cut protection and a soft and cool feel to the worker
Engineered to perform in combination with any cut resistant glove
11-280 Sleeves feature a seamless design and options of lengths and widths to ensure an excellent fit for all-day comfort
Available in three versions: Extra Wide, Wide, Narrow.

Technical Specification :

11-280 Sleeves



Laboratory Instrument India