



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
11-280 Sleeves

**Product Code :**  
LBNY-00017-11000522



**Description :**

INTERCEPT Cut Resistance Technology provides best-in-class cut protection and a soft and cool feel to the worker  
Engineered to perform in combination with any cut resistant glove  
11-280 Sleeves feature a seamless design and options of lengths and widths to ensure an excellent fit for all-day comfort  
Available in three versions: Extra Wide, Wide, Narrow.

**Technical Specification :**

11-280 Sleeves



**Laboratory Instrument India**