



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
School First Aid Kit

**Product Code :**  
LBNY-0003-1010011



**Description :**

School First Aid Kit

**Technical Specification :**

First aid kit box made from GI metal with red & ivory 7-layer powder coating which is rust proof, scratch-proof & corrosion free for long life.

Keep your first aid and medical supplies organized and easily accessible in this First Aid Kit and Medicine Storage Box.

Made with Front Transparent Acrylic Glass and Multi Compartment 3-Shelf for Better Organize of Medicines with Security Lock on Door for Safety

This First AID Kit Box is perfect for job sites, large office settings, or manufacturing facilities.

Suitable Homes, Offices, Hospitals, Schools, Police and Fire Departments, Construction Sites, Manufacturing Facilities, Sports Activities, Outdoor Activities, and Personal or Commercial Transportation.

Complete kit, plastic deluxe.

Size: 280 x 440 x 140mm mountable to the wall.



**Laboratory Instrument India**

