



Email : sales@laboratoryinstrumentindia.com

Product Name :
CPP Supported Black Neoprene Aprons

Product Code :
LBNY-00017-11000527



Description :

Provide excellent chemical resistance
All edges hemmed
Lightweight with four metal grommets, neck string, and waist ties
Top 2 in. (5cm) of bib is a double thickness of material
Dimensions: 45L x 35 in.W (114 x 89cm)
Thickness: 12.5 oz.

Technical Specification :

CPP Supported Black Neoprene Aprons



Laboratory Instrument India