



Email : sales@laboratoryinstrumentindia.com

Product Name :
G Clamp

Product Code :
LBNY-0003-1010012



Description :

G Clamp

Technical Specification :

The c-clamp is made of industrial-grade cast iron with professional heat treatment.

The c-shaped part is covered with red spaying plastic, and the T bar handle is chrome-plated. The well-made clamp is sturdy, rustproof, and durable.

ADJUSTABLE JAW OPENING SIZE: It is easy and quick to adjust the jaw opening size of the c-clamp set. Simply turn the T bar handle clockwise to narrow the jaw opening size and turn the opposite direction to enlarge the size.

T BAR HANDLE DESIGN: The heavy-duty c-clamp is designed with a T bar handle. There are threads on the T bar, which provides large torque and strong clamping force. The swivel pad on the tip increases the contact area and ensures that the clamping force is evenly distributed.

PRACTICAL TOOL: The 2-inch c-clamp is suitable for construction, metalworking, home improvement, general repairing and maintenance, woodworking, etc.

Plated steel with swivel pad screw maximum opening 75mm.



Laboratory Instrument India

