



Email : sales@laboratoryinstrumentindia.com

Product Name :

CPF 3 HD Personal-Protection Garments

Product Code :

LBNY-00017-11000533



Description :

Integrated gloves eliminate need for taping and provide improved protection
Heavy-duty front zipper closure and double storm flaps with wide hook and loop closure promote improved performance
Tactile glove integrated into garment†—Stretch Chloroprene with inner glove for wide spectrum chemical protection
Bigger zipper pull ensures easier donning and doffing
Garments deliver proven barrier protection against 110 toxic industrial chemicals and are designed to stand up to chemical warfare agents including Sarin, Mustard, and Lewisite † Choose from two styles: with three-panel respirator-fitted hood or with collar for compatibility with hooded PAPR.
Includes attached socks‡ with boot splash cover
Taped seams
Light tan.

Technical Specification :

CPF 3 HD Personal-Protection Garments



Laboratory Instrument India

