



Email : sales@laboratoryinstrumentindia.com

Product Name :
Plus 198 Coveralls

Product Code :
LBNY-00017-11000537



Description :

Plus garments are composed of flash spun high density polyethylene which creates a unique, nonwoven material that provides an ideal balance of protection, durability and comfort
Fabric and seams offer chemical permeation barrier to low concentration water-based inorganic chemicals
Serged and over-taped seams for protection and strength
Self-adhesive chin flap for tight seal of suit to the mask
Elastic face, wrists and ankles as well as glued-in waist elastic
Chemical protective clothing, Category III, Type 4-B, 5-B and 6-B
(protection against radioactive contamination)
Protection against infective agents including resistance to penetration by blood and body fluids using synthetic blood.
Coveralls provide an effective barrier against many water-based inorganic chemicals in low concentration and particles (down to 1.0 micron in size)Elastic thumb loops keep sleeves in place
White.

Technical Specification :

Plus 198 Coveralls

Laboratory Instrument India

