



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Plus 198 Coveralls

**Product Code :**  
LBNY-00017-11000537



**Description :**

Plus garments are composed of flash spun high density polyethylene which creates a unique, nonwoven material that provides an ideal balance of protection, durability and comfort  
Fabric and seams offer chemical permeation barrier to low concentration water-based inorganic chemicals  
Serged and over-taped seams for protection and strength  
Self-adhesive chin flap for tight seal of suit to the mask  
Elastic face, wrists and ankles as well as glued-in waist elastic  
Chemical protective clothing, Category III, Type 4-B, 5-B and 6-B  
(protection against radioactive contamination)  
Protection against infective agents including resistance to penetration by blood and body fluids using synthetic blood.  
Coveralls provide an effective barrier against many water-based inorganic chemicals in low concentration and particles (down to 1.0 micron in size)Elastic thumb loops keep sleeves in place  
White.

**Technical Specification :**

Plus 198 Coveralls

---

# Laboratory Instrument India

