



Email : sales@laboratoryinstrumentindia.com

Product Name :

600 Coverall with Hood, Elastic Ankles, Thumb loops

Product Code :

LBNY-00017-11000539



Description :

600 Coverall coveralls are ideal for occupations that require effective whole suit protection offered by a Type 4 suit because they are exposed to heavy liquid aerosols or toxic airborne solid particulates. 600 Coverall with Hood provides an effective barrier against low concentration water-based inorganic chemicals, small sized hazardous particles down to one micron, and biological hazards.

Enhanced design for increased flexibility and comfort:

New hood design: optimized to fit the contours of the face, ensuring fewer leakage points during head movement

Improved garment fit: extra flexibility and freedom of movement

Larger zipper puller: much easier to grasp when wearing gloves.

Technical Specification :

600 Coverall with Hood, Elastic Ankles, Thumb loops



Laboratory Instrument India