



Email : sales@laboratoryinstrumentindia.com

Product Name :
Isolation Gown

Product Code :
LBNY-00017-11000547



Description :

Designed to protect healthcare professionals and patients against cross-contamination
Features comfortable fit with thumb loop cuffs
Sewn seams provide for superior strength
These Isolation Gowns feature thumb loop cuffs to prevent sleeves from riding up and sewn seams to provide superior strength. Non-ISO Rated.
Single-use disposable
Made of spunbond polypropylene (latex-free).
Available in assorted sizes.

Technical Specification :

Isolation Gown



Laboratory Instrument India