



Email : sales@laboratoryinstrumentindia.com

Product Name :
Critical Cover Barrier Tech Gown

Product Code :
LBNY-00017-11000549



Description :

Durable and fluid impervious solution for lab and/or production areas
Provides protection for the wearer from particles, fluids and light chemical splash
Serged seams and elastic wrists.
Low particle count for use in controlled environments
Banded ties attached to the neck and waist for easy and secure closure.

Technical Specification :

Critical Cover Barrier Tech Gown



Laboratory Instrument India