



Email : sales@laboratoryinstrumentindia.com

Product Name :
Berkshire Choice Cleanroom Socks

Product Code :
LBNY-00017-11000552



Description :

99% Continuous filament polyester, 1% spandex
Mid-calf length with elastic cuff for secure fit
Soft and comfortable
One size fits all
Laundered and packaged in an ISO Class 5 cleanroom
Color: White
10 pairs per pack

Advantages

Comfortable and durable for all day wear
A cost effective, clean and disposable alternative to operators own socks
Autoclavable (121°C, 30 minutes)
Bulk packaging available.
Cleanroom laundry process reduces fiber and particle shedding and minimizes bioburden
Over 1000 times cleaner than regular polyester/cotton socks.

Technology

Cleanroom laundered and packaged
Constructed with continuous filament yarns to reduce fiber and particle shedding.

Technical Specification :

Berkshire Choice Cleanroom Socks



Laboratory Instrument India