



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Magnetic Compass Needle

**Product Code :**  
LBNY-0003-1010014



**Description :**

Magnetic Compass Needle

**Technical Specification :**

For studying magnetic phenomena such as operation of a compass, testing of magnetic poles, finding the direction of electric current etc.

Consists of a magnetised steel needle about 150mm long with brass bearing pivoted on a non-magnetic metal stand.

North pole is marked in red.

The needle when not in use is reversed and kept in main rod for safety.

On stand



**Laboratory Instrument India**