



Email : sales@laboratoryinstrumentindia.com

Product Name :
Magnetic Compass Needle

Product Code :
LBNY-0003-1010014



Description :

Magnetic Compass Needle

Technical Specification :

For studying magnetic phenomena such as operation of a compass, testing of magnetic poles, finding the direction of electric current etc.

Consists of a magnetised steel needle about 150mm long with brass bearing pivoted on a non-magnetic metal stand.

North pole is marked in red.

The needle when not in use is reversed and kept in main rod for safety.

On stand



Laboratory Instrument India