



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
BCR Cleanroom Socks

**Product Code :**  
LBNY-00017-11000558



**Description :**

Microdenier technology provides the sock with ideal cool feeling and comfort  
Provides microdenier softness for operator comfort  
Designed for all-day comfort  
Because of the enhanced wicking function, contact discomfort associated with wearing regular cotton or polyester socks is greatly reduced  
Reduced odor caused by foot perspiration  
Reduced sock degradation over extended wear due to oils and sweat  
Long lasting freshness and cooling for operator.

**Technical Specification :**

BCR Cleanroom Socks



**Laboratory Instrument India**