



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Berkshire BCR Cleanroom Socks

**Product Code :**  
LBNY-00017-11000559



**Description :**

Autoclavable (121°C, 30 minutes)  
96% microdenier nylon and 4% Lycra  
BCR Cleanroom socks are over 1,000 times cleaner than regular cotton/polyester socks  
One size fits all  
Extremely soft and comfortable  
White.

**Technical Specification :**

Berkshire BCR Cleanroom Socks



**Laboratory Instrument India**