




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Inter Metro PPE Rack

**Product Code :**  
LBNY-00017-11000564



**Description :**

PPE Rack features four-sided access to increase availability of sanitizer and PPE in high traffic areas. Shelves provide flexibility to station back-up supplies, other PPE, hold a waste container, or can be used for signage.

Exceptionally flexible

Can be reconfigured in one, two, or three-side configurations.

Sanitizer and PPE not included.

**Technical Specification :**

Inter Metro PPE Rack



**Laboratory Instrument India**