



Email : sales@laboratoryinstrumentindia.com

Product Name :
Heavyweight Terrycloth Heat Resistant Sleeves

Product Code :
LBNY-00017-11000565



Description :

Heavyweight Terrycloth Heat Resistant Sleeves with loop out terrycloth construction provide protection from temperatures up to 350°F.

Knit wrist helps prevent dirt and debris from entering the sleeve.

Protection from temperatures up to 350°C

Elastic top secures sleeve to upper arm.

Technical Specification :

Heavyweight Terrycloth Heat Resistant Sleeves



Laboratory Instrument India