



Email : sales@laboratoryinstrumentindia.com

Product Name :
Reusable Pot Sink Aprons

Product Code :
LBNY-00017-11000570



Description :

Suitable for food handling and preparation, general maintenance and chemical processing
Comfortable to wear with superior flexibility
Excellent resistance to punctures and abrasion, chemicals, fats, oils and grease
8 mil PVC/vinyl apron with neck ties
Individually packaged
White.

Technical Specification :

Reusable Pot Sink Aprons



Laboratory Instrument India