



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Reusable Pot Sink Aprons

**Product Code :**  
LBNY-00017-11000570



**Description :**

Suitable for food handling and preparation, general maintenance and chemical processing  
Comfortable to wear with superior flexibility  
Excellent resistance to punctures and abrasion, chemicals, fats, oils and grease  
8 mil PVC/vinyl apron with neck ties  
Individually packaged  
White.

**Technical Specification :**

Reusable Pot Sink Aprons



**Laboratory Instrument India**