



Email : sales@laboratoryinstrumentindia.com

Product Name :
PIP Neoprene Apron

Product Code :
LBNY-00017-11000571



Description :

Ideal for chemical plants, plating facilities and general maintenance
Provides superior protection from a wide range of chemicals including acids, caustics, solvents and oils
Lightweight black neoprene apron with tie strings
Thickness: 18-20 mil
Rustproof grommets
Reinforced bib
Hemmed overedge.

Technical Specification :

PIP Neoprene Apron



Laboratory Instrument India