




Email : sales@laboratoryinstrumentindia.com

Product Name :
SP Scienceware Vikem Vinyl Aprons

Product Code :
LBNY-00017-11000572



Description :

Provides protection to personnel in hospitals, schools, laboratories and industrial plants.
Strong tie backs double reinforced with metal grommets through three layers of material at waist
Made from 4mil thick vinyl
Stitched cotton binding securing all edges
Lightweight and comfortable
Non-absorbent and chemically resistant.

Technical Specification :

SP Scienceware Vikem Vinyl Aprons



Laboratory Instrument India