



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Cut Resistant Sleeve, Heavy Duty

**Product Code :**  
LBNY-00017-11000574



**Description :**

The white armguard is the highest cut resistance sleeve offered in the sleeve line. Constructed from a combination of high-performance materials offering A8 ANSI cut protection. Heavy Duty Cut Resistant Sleeve offers superior cut protection from arm lacerations. Prewashed and preshrunk for a close fit. Armguard offers great protection for the arm. Handguard II yarn technology. Outstanding cut and wear resistance for a long service life. Available in different lengths. 20 in. and 22 in. has a clip attached on the upper arm that allows the sleeve to stay in place.

**Technical Specification :**

Cut Resistant Sleeve, Heavy Duty



**Laboratory Instrument India**