



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cut Resistant Sleeve, Medium Duty

Product Code :
LBNY-00017-11000575



Description :

The medium duty armguard offers great protection for the arm and is available in different lengths. Gray armguard provides strength, and good cut resistance while hiding dirt
Features High Performance Fibers and Stainless Steel
Medium Duty Cut Resistant Sleeves offer superior cut protection from arm lacerations.
Constructed from a combination of high-performance materials offering A7 ANSI cut protection.
20 in. and 22 in. has a clip attached on the upper arm that allows the sleeve to stay in place.
Outstanding wear for a long service life
Prewashed and preshrunk for a close fit.

Technical Specification :

Cut Resistant Sleeve, Medium Duty



Laboratory Instrument India