



Email : sales@laboratoryinstrumentindia.com

Product Name :
Blue Denim Bib Aprons with Pockets

Product Code :
LBNY-00017-11000579



Description :

Ideal for industrial maintenance, construction, manufacturing and carpentry.
100% cotton construction and comfort
Metal grommets to hold neck strap and waist ties.
Serged seams and pre-shrunk for good launderability
Serged seams and pre-shrunk for good launderability
Available in a variety of pocket configurations.

Technical Specification :

Blue Denim Bib Aprons with Pockets



Laboratory Instrument India