



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
PIP Temp-Gard Extreme Temperature Apron

**Product Code :**  
LBNY-00017-11000584



**Description :**

Black silicone apron protects torso and upper legs in extreme temperature conditions  
Protects from splashes  
Easy to clean; will not absorb chemicals and will not stain  
Quick-release buckles at neck and waist for added safety  
Multi-layered liquid-proof black silicone fabric  
Impervious to cryogenic liquids and most other liquids and chemicals  
Does not become slippery, even at extreme temperatures  
Temperature range: -196° to +260°C (-321° to +500°F)  
Cryogenic cleaning, cryogenic systems maintenance, conservation of food produce.  
Superior flexibility; will not stiffen.

**Technical Specification :**

PIP Temp-Gard Extreme Temperature Apron



**Laboratory Instrument India**

