



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Protective Products ATA Cut Resistant Sleeve 1-ply Design

**Product Code :**  
LBNY-00017-11000590



**Description :**

Cut Resistant Sleeve delivers ANSI A4 cut protection in a remarkably thin, 1-ply design. Designed for outstanding stretch and fabric memory, a 2.5 sleeve fits just about any size arm, big or small  
Lightweight 1-Ply Seamless Knit  
With stretch and fabric memory, it fits just about any size arm, and can be worn over clothing or against the skin.  
18" length includes a Thumb Hole to hold the sleeve in place, keep the wrist protected  
Machine washable.  
Get solid A4 cut protection in a remarkably thin, light weight sleeve  
Cool and comfortable on the arm; it can be worn over clothing or directly against the skin.

**Technical Specification :**

Protective Products ATA Cut Resistant Sleeve 1-ply Design



**Laboratory Instrument India**