



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Micrometer Screw Gauge

**Product Code :**  
LBNY-0003-1010018



**Description :**

Micrometer Screw Gauge

**Technical Specification :**

Made of nickel plated brass, with ratchet top, accurately machined stainless steel rod.

Supplied in case.

Accuracy is perfect for a hobbyist or student.

Range: 0-25 x 0.01mm.

Supplied with a protective case.

Ratchet of iron/steel.



**Laboratory Instrument India**