



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
PIP Clear Polyethylene Apron

**Product Code :**  
LBNY-00017-11000596



**Description :**

Ideal for food handling and preparation and general maintenance applications.  
Economical protection for light- to medium-duty applications  
Excellent resistance to chemicals, fats, oils and grease  
Lightweight embossed 20mil polyethylene apron with ties at neck  
Individually packaged  
Easy to clean  
Die cut for consistent sizing.

**Technical Specification :**

PIP Clear Polyethylene Apron



**Laboratory Instrument India**