



Email : sales@laboratoryinstrumentindia.com

Product Name :
PIP Clear Polyethylene Apron

Product Code :
LBNY-00017-11000596



Description :

Ideal for food handling and preparation and general maintenance applications.
Economical protection for light- to medium-duty applications
Excellent resistance to chemicals, fats, oils and grease
Lightweight embossed 20mil polyethylene apron with ties at neck
Individually packaged
Easy to clean
Die cut for consistent sizing.

Technical Specification :

PIP Clear Polyethylene Apron



Laboratory Instrument India