



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
PIP Hi-Vis Vests

**Product Code :**  
LBNY-00017-11000604



**Description :**

Side elastic for improved comfort and fit

Breathable polyester mesh fabric

Hook and loop closure

Ideal for additional levels of safety for low risks activities where ANSI compliance is not required, such as shopping cart retrievers, dog walkers, parking lot attendants, joggers, delivery vehicle drivers, recreational workers.

Available in hi-vis lime yellow and hi-vis orange.

**Technical Specification :**

PIP Hi-Vis Vests



**Laboratory Instrument India**