



Email : sales@laboratoryinstrumentindia.com

Product Name :
Uniform Technology Coverall Grid, White, Long

Product Code :
LBNY-00017-11000606



Description :

Ideal for General CE Spaces/Applications (Non-ESD sensitive Medical Device, Electronics & Aerospace, Food Processing/Handling or Packaging, Cultivation Industry, University Labs, Imaging/Film coating Manufactures, General Medical/Dental Facilities)

Long style cleanroom coveralls feature a military collar, raglan sleeves, knit cuffs, and a coil zipper with top snaps.

Coverall

Knit Cuffs

Coil Zipper with Top Snaps.

Military Collar

Raglan Sleeves.

Technical Specification :

Uniform Technology Coverall Grid, White, Long



Laboratory Instrument India

